

# TRAUFIX™




**TRAUFIX**  
READY TO RUN

CANULATED SUBTALAR SCREW FOR CORRECTION OF FLAT FOOT VALGUS

# TARFIX CANULATED SUBTALAR SCREW FOR CORRECTION OF FLAT FOOT VALGUS


## CONICAL DESIGN

- ▶ Press-fit; stable hold.
- ▶ Anatomical adaptation.
- ▶ Easy insertion and removal.




## 4 ROUNDED CHANNELS

- ▶ Reduced direct contact area: less friction and less irritation.
- ▶ Promotes tissue integration and blood irrigation.
- ▶ Controlled structural flexibility.



## TWO HOLES IN THE BODY

improved osseointegration



## INTERNAL THREADS IN THE SCREW BODY

Allows the extractor screwdriver (251.01) to be inserted.  
Opposite thread for efficient extraction.



## MEASURES



8 mm

9 mm

10 mm

11 mm

12 mm


13 mm

14 mm



CROSS-SHAPED SCREWDRIVER ENTRY  
to improve traction during insertion.

### Specifications



### Length Table

TARFIX CANULATED SUBTALAR SCREW FOR CORRECTION OF FLAT FOOT VALGUS

CODE	DIAMETER	LT (mm)
<b>262.08</b>	8	14.00
<b>262.09</b>	9	14.00
<b>262.10</b>	10	14.00
<b>262.11</b>	11	16.00
<b>262.12</b>	12	16.00
<b>262.13</b>	13	19.00
<b>262.14</b>	14	19.00



Titanium Alloy 6Al4V ELI ASTM F136

Special for surgical implants

Design by **TRAUFIX**  
READY TO RUN



MADE IN MEXICO

# TARFIX SET

## QTY. EQUIPMENT

- 1 Tarfix screw removal tool. (Cód. 251.01)
- 1 Crosshead screwdriver for inserting Tarfix screws. (Cód. 251.02)
- 1 Blunt-tipped guide wire, diameter 2.0mm, length 265mm. (Cód. 318.2026)
- 1 Tarfix screw tester, sizes:
  - 8mm (Code 319.08), 9mm (Code 319.09), 10mm (Code 319.10), 11mm (Code 319.11),
  - 12mm (Code 319.12), 13mm (Code 319.13) and 14mm (Code 319.14).



Our containers for implants and screws are manufactured with **TECAPRO MT**  
(Heat Stabilized Polypropylene) PP Medical Grade

Ideal material for sterilization containers

### TRAUFIX

internationalsales@fixier.com.mx  
(+52) 55 8106 3240

[www.traufix.com](http://www.traufix.com)

### DIPROMEDIC MEXICO

customerservice@dipromedic.com  
(55) 9000 4888/ 89/ 90/ 91

### DIPROMEDIC PERU

gustavo.dominguez@dipromedic.com  
ventas.interdipro.peru@dipromedic.com  
(+51) 947 732 069

