


Intralock


Titanium Tibia




TRAUFIX
READY TO RUN

INTRALOCK TIBIAL NAIL


Two options of mid-lateral locking allows primary compression or controlled secondary dynamization.




DISTAL HOLE FOR DEVICE




Innovative distal locking device for precise collocation of 4.5mm screw



TO PERFORM LOCKING WITH DEVICE
PRELIMINARY STEP. Calibration of distal strip
to generate hole of distal locking device.





Distal and proximal locking system with multiple options, prevent axial rotation, improve stability and ensure fracture reduction.



Code 166



Code 169.11

MEASURES

TIBIAL INTRAMEDULLARY NAIL INTRALOCK

Length	Diameter 8.5mm	Diameter 9.5mm
255mm	Code 163.255	Code 162.255
270mm	Code 163.270	Code 162.270
285mm	Code 163.285	Code 162.285
300mm	Code 163.300	Code 162.300
315mm	Code 163.315	Code 162.315
330mm	Code 163.330	Code 162.330
345mm	Code 163.345	Code 162.345
360mm	Code 163.360	Code 162.360
375mm	Code 163.375	Code 162.375

4.5mm LOCKING SCREW FOR INTRAMEDULLARY NAIL INTRALOCK (CODE 166, LENGTH 25mm-80mm)
 END CAP FOR TIBIAL NAIL INTRALOCK (CODE 169.11)

Titanium Alloy 6AL4V ELI ASTM F136

Special for surgical implants

Design by **TRAUFIX**
 READY TO RUN

MADE IN MEXICO






TIBIAL NAIL SET

INSTRUMENTATION AND IMPLANTS SET. ALL-IN-ONE

QTY. EQUIPMENT

- 1 T-handle for quick coupling
- 1 Cannulated Starter Drill Ø12.5
- 1 T-handle with hex blade Ø6.5
- 1 Initiator punch
- 2 Nail Extractor Connector
- 1 Compression screw
- 1 Drill Ø5.2
- 2 Graduated Drill Ø3.2
- 1 Graduated drill with cap Ø3.2
- 3 Threaded tip guide Ø3.0
- 2 Drill Ø4.0
- 1 Drill with cap Ø4.0
- 2 Flexible rhyme
- 2 Nail Fastener Screw
- 1 Fabric protector
- 1 Drill Guide Ø 12
- 1 Spanish key Ø 11
- 1 Lock tester
- 1 Big Allen Key Ø 5
- 1 Small Allen Key Ø 3
- 1 Rule
- 1 Key Ø4
- 1 Drill sleeve interior and exterior Ø 4 y Ø 8
- 1 Drill sleeve interior and exterior Ø 4 y Ø 8
- 1 Drill sleeve exterior Ø 5.2
- 1 Tissue Dilator(Punch) Ø 8.0
- 1 Thin rope distal rail screw Ø 8.5 Ø 3.2 and Ø 3.2
- 1 Depth meter 90
- 1 Stylus with T-handle
- 1 T-handle drill bit Ø 5
- 1 Guide Wire Guard
- 1 Nail Holder for Extraction SW 3.5
- 1 Hex T-Handle Bolt Driver Ø 4.5
- 1 T-Handle Hex Driver Ø 6.5
- 1 AO Hex Bolt Driver Bit 4.5 SW3.5
- 1 Proximal locking device
- 1 Distal locking device
- 1 Distal Locking Device for Stylus
- 1 Nail Extraction Guide
- 2 Locating screw
- 1 Screw strip
- 1 Nail holder for extraction
- 1 Rhyming Tips Set
Ø 7.5 Ø 8.0 Ø 8.5 Ø 9.0 Ø 9.5 Ø 10 Ø 10.5
Ø 11 Ø 11.5 Ø 12
- 1 Guide Wire Handle
- 1 Blocking strip U shape
- 1 Insertion handle
- 1 Cannulated chisel



TRAUFIX

internationalsales@fixier.com.mx
(+52) 55 8106 3240

www.traufix.com

DIPROMEDIC MEXICO

customerservice@dipromedic.com
(55) 9000 4888/ 89/ 90/ 91

DIPROMEDIC PERU

gustavo.dominguez@dipromedic.com
ventas.interdipro.peru@dipromedic.com
(+51) 947 732 069

 traufix.mx

 traufix